

OBTURATOR VISUAL FOR RESECTOSCOPY

27 / 28.5Fr.; with 7Fr. working channel

SKU: 630-230-27

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Obturators

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Obturators