

OCHSNER AORTIC CLAMP

angled; atraumatic

22.0cm (8.6")



SKU: 142-010-23

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Atrauma Forceps
Aesculap	FB495
KLS Martin	24-317-23
AUTHORS	OCHSNER
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Atrauma Forceps
Aesculap	FB495
KLS Martin	24-317-23
AUTHORS	OCHSNER

SPECIFICATIONS

[angled](#)