

## O´CONNOR HOOK

sharp

16.0cm (6.4")

**SKU:** 062-074-16

**Categories:** [Open Surgery](#), [Retractors](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>KLS Martin</b>	15-070-16
<b>AUTHORS</b>	<a href="#">O´CONNOR</a>
<b>SPECIFICATIONS</b>	<a href="#">sharp</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>KLS Martin</b>	15-070-16
<b>AUTHORS</b>	<a href="#">O´CONNOR</a>
<b>SPECIFICATIONS</b>	<a href="#">sharp</a>