

OLIVECRONA TRIGEMINAL SCISSORS

serrated

working length: 120 mm




SKU: 010-469-23

Categories: [Open Surgery](#), [Scissors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Scissors for Neuro Surgery , Scissors for Vascular Surgery
Aesculap	BC665
KLS Martin	11-475-23
AUTHORS	OLIVECRONA
Medicon	57.25.70
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Scissors for Neuro Surgery , Scissors for Vascular Surgery
Aesculap	BC665
KLS Martin	11-475-23

AUTHORS

[OLIVECRONA](#)

Medicon

57.25.70

SPECIFICATIONS

[straight](#)