

OLSEN-HEGAR NEEDLE HOLDER

straight; serrated; with scissors; groove

19.0cm (7.6")


SKU: 050-310-19

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM291
KLS Martin	20-464-18
AUTHORS	HEGAR , OLSEN
Medicon	10.60.19
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM291
KLS Martin	20-464-18

AUTHORS

[HEGAR, OLSEN](#)

Medicon

10.60.19

SPECIFICATIONS

[serrated](#), [straight](#)