

OLSEN-HEGAR NEEDLE HOLDER

straight; serrated; with scissors; groov

14.0cm (5.6")


SKU: 050-310-14

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM290
KLS Martin	20-464-14
AUTHORS	HEGAR , OLSEN
Medicon	10.60.14
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM290
KLS Martin	20-464-14

AUTHORS

[HEGAR](#), [OLSEN](#)

Medicon

10.60.14

SPECIFICATIONS

[serrated](#), [straight](#)