

## OLSEN-HEGAR NEEDLE HOLDER

straight; serrated; with scissors; groove

14.0cm (5.6")


**SKU:** 050-310-14

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                |
| <b>TYPE</b>           | <a href="#">Needle Holders</a>                      |
| <b>Aesculap</b>       | BM290   |
| <b>KLS Martin</b>     | 20-464-14   |
| <b>AUTHORS</b>        | <a href="#">HEGAR</a> , <a href="#">OLSEN</a>       |
| <b>Medicon</b>        | 10.60.14  |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                   |                                |
|-------------------|--------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>           |
| <b>TYPE</b>       | <a href="#">Needle Holders</a> |
| <b>Aesculap</b>   | BM290                          |
| <b>KLS Martin</b> | 20-464-14                      |

**AUTHORS**

[HEGAR](#), [OLSEN](#)

**Medicon**

10.60.14

**SPECIFICATIONS**

[serrated](#), [straight](#)