


## ORINGER LIP RETRACTOR

fig. 3

11.0cm (4.4")



**SKU:** 114-243-03

**Categories:** [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)

### ZUSÄTZLICHE INFORMATIONEN

|                   |   |
|-------------------|---|
| <b>UNIT</b>       | <a href="#">pcs.</a>  |
| <b>TYPE</b>       | <a href="#">Cheek and Lip Retractors, Tongue Depressors</a> |
| <b>KLS Martin</b> | 38-240-03   |
| <b>AUTHORS</b>    | <a href="#">ORINGER</a>                                     |

### ZUSÄTZLICHE INFORMATIONEN

|                   |   |
|-------------------|---|
| <b>UNIT</b>       | <a href="#">pcs.</a>  |
| <b>TYPE</b>       | <a href="#">Cheek and Lip Retractors, Tongue Depressors</a> |
| <b>KLS Martin</b> | 38-240-03   |
| <b>AUTHORS</b>    | <a href="#">ORINGER</a>                                     |