


ORINGER LIP RETRACTOR

for children

9.5cm (3.8")



SKU: 114-243-00

Categories: [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Cheek and Lip Retractors, Tongue Depressors
KLS Martin	38-240-00
AUTHORS	ORINGER

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Cheek and Lip Retractors, Tongue Depressors
KLS Martin	38-240-00
AUTHORS	ORINGER