

## OTOLOGY SUCTION TUBE

tapered; curved; Luer connector

Ø 3.0mm; WL: 17.5cm

**SKU:** 086-504-30

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



## ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Suction Tubes](#)

## ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Suction Tubes](#)