

## OTOLOGY SUCTION TUBE

cylindrical; curved; Luer connector

Ø 5.0mm; WL: 11.0cm

**SKU:** 086-511-50

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>