

OUTER SHEATH; 18.5FR.; ALLROUND PERF

for inner sheath 582-180-26; rotatable

QL connection, colour code green



SKU: 582-180-24

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Mini Resectoscopy (Midi)
TYPE	Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Mini Resectoscopy (Midi)
TYPE	Sheaths