

## OUTER SHEATH; 18.5FR.; ALLROUND PERF

for inner sheath 582-180-26; rotatable

QL connection, colour code green



**SKU:** 582-180-24

**Categories:** [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Mini Resectoscopy (Midi)</a>
<b>TYPE</b>	<a href="#">Sheaths</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Mini Resectoscopy (Midi)</a>
<b>TYPE</b>	<a href="#">Sheaths</a>