

## OVERHOLT-FINO LIGATURE FORCEPS

curved; serrated

23.0cm (9.2")




**SKU:** 040-426-03

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dissecting and Ligature Forceps</a>
<b>KLS Martin</b>	13-448-03
<b>AUTHORS</b>	<a href="#">FINO</a> , <a href="#">OVERHOLT</a>
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dissecting and Ligature Forceps</a>
<b>KLS Martin</b>	13-448-03
<b>AUTHORS</b>	<a href="#">FINO</a> , <a href="#">OVERHOLT</a>
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>