

## OVERHOLT-FINO LIGATURE FORCEPS

curved; serrated

21.5cm (8.6")




**SKU:** 040-426-00

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                              |
| <b>TYPE</b>           | <a href="#">Dissecting and Ligature Forceps</a>   |
| <b>Aesculap</b>       | BJ080   |
| <b>KLS Martin</b>     | 13-448-00   |
| <b>AUTHORS</b>        | <a href="#">FINO</a> , <a href="#">OVERHOLT</a>   |
| <b>SPECIFICATIONS</b> | <a href="#">curved</a> , <a href="#">serrated</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                   |   |
|-------------------|---|
| <b>UNIT</b>       | <a href="#">pcs.</a>                            |
| <b>TYPE</b>       | <a href="#">Dissecting and Ligature Forceps</a> |
| <b>Aesculap</b>   | BJ080   |
| <b>KLS Martin</b> | 13-448-00                                       |
| <b>AUTHORS</b>    | <a href="#">FINO</a> , <a href="#">OVERHOLT</a> |

## SPECIFICATIONS

[curved](#), [serrated](#)