

OVERHOLT LIGATURE FORCEPS

curved; serrated

29.5cm (11.5")




SKU: 040-427-00

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
Aesculap	BJ090
KLS Martin	13-449-00
AUTHORS	OVERHOLT
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
Aesculap	BJ090
KLS Martin	13-449-00
AUTHORS	OVERHOLT

SPECIFICATIONS

[curved](#), [serrated](#)