

OVERHOLT LIGATURE FORCPES

s-shaped; serrated

21.0cm (8.4")




SKU: 040-423-03

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
Aesculap	BJ017
KLS Martin	13-417-03
AUTHORS	OVERHOLT
Medicon	16.64.21
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
Aesculap	BJ017
KLS Martin	13-417-03

AUTHORS

[OVERHOLT](#)

Medicon

16.64.21

SPECIFICATIONS

[curved](#), [serrated](#)