

OVERHOLT LIGATURE FORCPES

s-shaped; serrated

21.5cm (8.6")




SKU: 040-423-04

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Dissecting and Ligature Forceps |
| Aesculap | BJ018 |
| KLS Martin | 13-417-04 |
| AUTHORS | OVERHOLT |
| Medicon | 16.64.22 |
| SPECIFICATIONS | curved , serrated |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---|
| UNIT | pcs. |
| TYPE | Dissecting and Ligature Forceps |
| Aesculap | BJ018 |
| KLS Martin | 13-417-04 |

AUTHORS

[OVERHOLT](#)

Medicon

16.64.22

SPECIFICATIONS

[curved](#), [serrated](#)