

## OVERHOLT LIGATURE FORCPES

s-shaped; serrated

21.5cm (8.6")


**SKU:** 040-423-04

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BJ018
<b>KLS Martin</b>	13-417-04
<b>AUTHORS</b>	<a href="#">OVERHOLT</a>
<b>Medicon</b>	16.64.22
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BJ018
<b>KLS Martin</b>	13-417-04

## AUTHORS

[OVERHOLT](#)

## Medicon

16.64.22

## SPECIFICATIONS

[curved](#), [serrated](#)