

OVERHOLT LIGATURE FORCPES

s-shaped; serrated

24.0cm (9.6")




SKU: 040-423-05

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
KLS Martin	13-417-05
AUTHORS	OVERHOLT
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
KLS Martin	13-417-05
AUTHORS	OVERHOLT
SPECIFICATIONS	curved , serrated