

OVERHOLT LIGATURE FORCPES

s-shaped; serrated

18.5cm (7.4")




SKU: 040-423-01

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
Aesculap	BJ015
KLS Martin	13-417-01
AUTHORS	OVERHOLT
Medicon	16.64.18
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
Aesculap	BJ015
KLS Martin	13-417-01

AUTHORS

[OVERHOLT](#)

Medicon

16.64.18

SPECIFICATIONS

[curved](#), [serrated](#)