

PEAN-NIPPON HAEMOSTATIC FORCEPS

straight; serrated

18.0cm (7.2")



SKU: 040-404-18

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
KLS Martin	13-380-18
AUTHORS	NIPPON , PEAN
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
KLS Martin	13-380-18
AUTHORS	NIPPON , PEAN
SPECIFICATIONS	serrated , straight