

## PEAN-NIPPON HAEMOSTATIC FORCEPS

straight; serrated

14.5cm (5.6")



**SKU:** 040-404-14

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>KLS Martin</b>	13-380-14
<b>AUTHORS</b>	<a href="#">NIPPON</a> , <a href="#">PEAN</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>KLS Martin</b>	13-380-14
<b>AUTHORS</b>	<a href="#">NIPPON</a> , <a href="#">PEAN</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>