

## PEAN-NIPPON HAEMOSTATIC FORCEPS

straight; serrated

14.5cm (5.6")



**SKU:** 040-404-14

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                |
| <b>TYPE</b>           | <a href="#">Artery Forceps</a>                      |
| <b>KLS Martin</b>     | 13-380-14   |
| <b>AUTHORS</b>        | <a href="#">NIPPON</a> , <a href="#">PEAN</a>       |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                |
| <b>TYPE</b>           | <a href="#">Artery Forceps</a>                      |
| <b>KLS Martin</b>     | 13-380-14   |
| <b>AUTHORS</b>        | <a href="#">NIPPON</a> , <a href="#">PEAN</a>       |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a> |