

PENNINGTON HAEMOSTATIC FORCEPS

straight; serrated

20.0cm (8")




SKU: 040-358-20

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	EA037
KLS Martin	13-280-20
AUTHORS	PENNINGTON
Medicon	17.34.20
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	EA037
KLS Martin	13-280-20

AUTHORS

[PENNINGTON](#)

Medicon

17.34.20

SPECIFICATIONS

[serrated](#), [straight](#)