

PENNINGTON HAEMOSTATIC FORCEPS

straight; serrated

15.5cm (6.2")

SKU: 040-358-15

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Artery Forceps](#)


Aesculap BH902

KLS Martin 13-280-15

AUTHORS [PENNINGTON](#)

Medicon 17.34.15

SPECIFICATIONS [serrated](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Artery Forceps](#)

Aesculap BH902

KLS Martin 13-280-15

AUTHORS

[PENNINGTON](#)

Medicon

17.34.15

SPECIFICATIONS

[serrated](#), [straight](#)