

PENNINGTON HAEMOSTATIC FORCEPS

straight; serrated

15.5cm (6.2")




SKU: 040-358-15

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH902
KLS Martin	13-280-15
AUTHORS	PENNINGTON
Medicon	17.34.15
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH902
KLS Martin	13-280-15

AUTHORS

[PENNINGTON](#)

Medicon

17.34.15

SPECIFICATIONS

[serrated](#), [straight](#)