

PLESTER SELFRETAINING RETRACTOR

right

11.0cm (4.4")



SKU: 070-673-11

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV091
KLS Martin	15-759-11
AUTHORS	PLESTER
Medicon	24.26.01



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV091
KLS Martin	15-759-11
AUTHORS	PLESTER

Medicon

24.26.01