

PLESTER SELFRETAINING RETRACTOR

sharp; 2X2 teeth

13.0cm (5.2")

SKU: 070-830-13

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Selfretaining Retractors |
| Aesculap | BV078 |
| KLS Martin | 15-756-13 |
| AUTHORS | PLESTER |
| SPECIFICATIONS | sharp , toothed |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Selfretaining Retractors |
| Aesculap | BV078 |
| KLS Martin | 15-756-13 |
| AUTHORS | PLESTER |

SPECIFICATIONS

[sharp, toothed](#)