

PLESTER SELFRETAINING RETRACTOR

sharp; 2X2 teeth

13.0cm (5.2")




SKU: 070-830-13

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV078
KLS Martin	15-756-13
AUTHORS	PLESTER
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV078
KLS Martin	15-756-13
AUTHORS	PLESTER

SPECIFICATIONS

[sharp, toothed](#)