

PLESTER SELFRETAINING RETRACTOR

blunt; 2X2 teeth

13.0cm (5.2")




SKU: 070-831-13

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV079
KLS Martin	15-755-13
AUTHORS	PLESTER
Medicon	24.26.03
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV079
KLS Martin	15-755-13

AUTHORS

[PLESTER](#)

Medicon

24.26.03

SPECIFICATIONS

[blunt, toothed](#)