

## PLESTER SELFRETAINING RETRACTOR

blunt; 2X2 teeth

13.0cm (5.2")




**SKU:** 070-831-13

**Categories:** [Open Surgery](#), [Retractors](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>Aesculap</b>	BV079
<b>KLS Martin</b>	15-755-13
<b>AUTHORS</b>	<a href="#">PLESTER</a>
<b>Medicon</b>	24.26.03
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>Aesculap</b>	BV079
<b>KLS Martin</b>	15-755-13

## AUTHORS

[PLESTER](#)

Medicon

24.26.03

## SPECIFICATIONS

[blunt, toothed](#)