

PLESTER SUCTION TUBE

6 Fr.; 30° angled

14.0cm (5.2") working length

SKU: 085-032-06

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
Aesculap	GF940
AUTHORS	PLESTER
Medicon	09.84.06

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
Aesculap	GF940
AUTHORS	PLESTER
Medicon	09.84.06