

PLESTER SUCTION TUBE

6 Fr.; 30° angled

14.0cm (5.2") working length

SKU: 085-032-06

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| Aesculap | GF940 |
| AUTHORS | PLESTER |
| Medicon | 09.84.06 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| Aesculap | GF940 |
| AUTHORS | PLESTER |
| Medicon | 09.84.06 |