

## PLESTER SUCTION TUBE

8 Fr.; 30° angled

14.0cm (5.2") working length

**SKU:** 085-032-08

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>Aesculap</b>	GF941
<b>AUTHORS</b>	<a href="#">PLESTER</a>
<b>Medicon</b>	09.84.08

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>Aesculap</b>	GF941
<b>AUTHORS</b>	<a href="#">PLESTER</a>
<b>Medicon</b>	09.84.08