


POLYPECTOMY LOOP

oval; 20mm; monofil

tubus (5fr.) 1.8 x 250cm; single use



SKU: PK-4EOPK1820250M-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)