


## POLYPECTOMY LOOP

oval; 20mm; monofil

tubus (5fr.) 1.8 x 250cm; single use



**SKU:** PK-4EOPK1820250M-S

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)