

POLYPECTOMY LOOP, MONOFIL, ROTATABLE

oval 30mm, insulated tip

2.6mm, length: 230cm



SKU: PK-4POPK2630230IMDZ

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Polypectomy Loops, EMR-Snares](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Polypectomy Loops, EMR-Snares](#)