


POLYPECTOMY LOOP , MONOFIL (SINGLE USE)

snare: oval; 30mm, box 5/pcs

tubus 2.3 x 230cm



SKU: PK-4EOPK2330230M-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Single Use Instruments](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Single Use Instruments](#)