


## POLYPECTOMY LOOP (SINGLE USE)

snare: hexagonal; 30mm, box 5/pcs

tubus 2.3 x 230cm



**SKU:** PK-4EHPK2330230-S

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)