


## POLYPECTOMY LOOP, SPRING FLAT REINFORCED

oval; 20mm, with teeth

tubus (7fr.) 2.3 x 250cm; single use



**SKU:** PK-4EOPK2330250FY-S

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)