


## POLYPECTOMY LOOP

hexagonal; 25mm

tubus (5fr.) 1.8x250cm; single use



**SKU:** PK-4EHPK1825250-S

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)