

## POOLE SUCTION TUBE; STRAIGHT

perforated and removable external tube

inner Ø 4.0mm; outer Ø 8.0mm; 24.5cm

**SKU:** 085-020-08

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>AUTHORS</b>	<a href="#">POOLE</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>AUTHORS</b>	<a href="#">POOLE</a>