

POWERGRIP/ORBITARIS ELECTRODE INSERT

34.0cm; grasping forceps with 2x3 teeth

with P-COAT-insulation; silver gray

SKU: 628-352-53

Categories: [Endoscopy](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Instruments, bipolar
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Instruments, bipolar
TYPE	Electrodes