

PRESBYTERIAN TUBING CLAMP

straight; smooth

18.0cm (7.2")

SKU: 032-012-18

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

AUTHORS [PRESBYTERIAN](#)

Medicon 17.72.18

SPECIFICATIONS [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

AUTHORS [PRESBYTERIAN](#)

Medicon 17.72.18

SPECIFICATIONS [straight](#)