

PRESBYTERIAN TUBING CLAMP

straight; smooth

18.0cm (7.2")




SKU: 032-012-18

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
AUTHORS	PRESBYTERIAN
Medicon	17.72.18
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
AUTHORS	PRESBYTERIAN
Medicon	17.72.18
SPECIFICATIONS	straight