

PRESBYTERIAN TUBING CLAMP

straight; smooth

20.0cm (8")




SKU: 032-012-20

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
AUTHORS	PRESBYTERIAN
Medicon	17.72.20
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
AUTHORS	PRESBYTERIAN
Medicon	17.72.20
SPECIFICATIONS	straight