

PRESBYTERIAN TUBING CLAMP

straight; smooth

20.0cm (8")

SKU: 032-012-20

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

AUTHORS [PRESBYTERIAN](#)

Medicon 17.72.20

SPECIFICATIONS [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

AUTHORS [PRESBYTERIAN](#)

Medicon 17.72.20

SPECIFICATIONS [straight](#)