


PROBE

Ø 2.0mm

16.0cm (6.4")



SKU: 080-004-16

Categories: [Open Surgery](#), [Probes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN136
KLS Martin	16-100-16
Medicon	09.05.16
SPECIFICATIONS	without eye

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN136
KLS Martin	16-100-16
Medicon	09.05.16

SPECIFICATIONS

[without eye](#)