


## PROBE

Ø 2.0mm

18.0cm (7.2")



**SKU:** 080-004-18

**Categories:** [Open Surgery](#), [Probes](#)

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Buttoned Probes</a>
<b>Aesculap</b>	BN138
<b>KLS Martin</b>	16-100-18
<b>Medicon</b>	09.05.18
<b>SPECIFICATIONS</b>	<a href="#">without eye</a>

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Buttoned Probes</a>
<b>Aesculap</b>	BN138
<b>KLS Martin</b>	16-100-18
<b>Medicon</b>	09.05.18

## SPECIFICATIONS

[without eye](#)